

# ROOT ELEVATORS



## VALUED BUYERS

Our Elevators feature all stainless steel construction for greater durability. This allows us to argon weld our tips to the handle-not simply press fit them as is done in brass handle elevators. This welding process gives us unsurpassed strength at the critical tip handle juncture. Each perfectly forged tip undergoes a heat treatment in our high-tech, computerized vacuum furnace to ensure optimal hardness without risk of brittleness. We regularly manufacturing below articles in showing handles & pattern A,B,C,D,E & F for better grip & perfection. We are also fully equipped & capable to develop these instruments as per your own specific requirements & demands in Materials, Sizes, Shapes and Patterns.



A



B



C



D



E



F



ROOT ELEVATORS 167 TO 186

ROOT TIP PICKS, SYNDESMOTOMES  
& PERIOTOMES 187 TO 189



## A.N.A. AMITY IMPEX

A.N.A Square, Airport Road Muradpur, Sialkot-51360, Pakistan.  
Tel: +92-52-4268654, 4564127 Fax: +92-52-4261027, Mob: +92-300-8618654.  
E-mai: info@anaamity.com, anaamityimpex@gmail.com, sales@anaamity.com  
Web : www.anaamity.com



ROOT ELEVATORS

*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*



**ANA-003-100**  
# 1  
straight

2.5 mm

3 mm



**ANA-003-101**  
Fig. 1B  
Bein

4 mm



**ANA-003-102**  
Fig. 2B  
Bein

3 mm



**ANA-003-103**  
Fig. 3B  
Bein

4 mm



**ANA-003-104**  
Fig. 4B  
Bein

3.2 mm



**ANA-003-105**  
Fig. 5B  
Bein

3.8 mm



**ANA-003-106**  
Fig. 6B  
Bein

5.0 mm



**ANA-003-107**  
Fig. 4  
straight

3 mm



**ANA-003-108**  
Fig. 1F  
Flohr

3.5 mm



**ANA-003-109**  
Fig. 2F  
Flohr



**ANA-003-110**  
Fig. 3F  
Flohr



**ANA-003-111**  
Fig. 4B SA  
Schmeckebier  
Apexo  
Concave tip



**ANA-003-112**  
Fig. 5B SA  
Schmeckebier  
Apexo  
Concave tip



**ANA-003-113**  
Fig. 4 SA  
Schmeckebier  
Apexo  
Flat tip




**ANA-003-114**  
Fig. 5 SA  
Schmeckebier  
Apexo  
Flat tip

# ROOT ELEVATORS



*Елеватори екстракційні  
Wurzelheber*

*Elévateurs à racines  
Elevadores para raíces*



ROOT ELEVATORS

*Елеватори екстракційні*  
*Wurzelheber*

*Élévateurs à racines*  
*Elevadores para raíces*



**ANA-003-123**  
Fig. 1K  
Kopp

3.6 mm



**ANA-003-124**  
Fig. 2K  
Kopp

4.2 mm



**ANA-003-125**  
Fig. 3K  
Kopp

4.2 mm



**ANA-003-126**  
Fig. 4K  
Kopp

3.8 mm



**ANA-003-127**  
Fig. 5K  
Kopp



**ANA-003-128**  
Fig. 6K  
Kopp

2.5 mm



**ANA-003-129**  
Fig. 1

2.5 mm



**ANA-003-130**  
Fig. 2

4.5 mm



**ANA-003-131**  
Fig. 3

4.5 mm



**ANA-003-132**  
Fig. 4

3.3 mm



**ANA-003-133**  
Fig. 1



**ANA-003-134**  
Fig. 2

# ROOT ELEVATORS



*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*



**ANA-003-135**  
Fig. 71  
Miller-Apexo



**ANA-003-136**  
Fig. 72  
Miller-Apexo



**ANA-003-137**  
Fig. 73  
Miller-Apexo



**ANA-003-138**  
Fig. 74  
Miller-Apexo



**ANA-003-139**  
Fig. 19



**ANA-003-140**  
Fig. 20



**ANA-003-141**  
Fig. 34S  
Seldin



**ANA-003-142**  
Fig. 34  
Seldin



**ANA-003-143**  
Fig. 2S  
Seldin



**ANA-003-144**  
Fig. 3S  
Seldin

ROOT ELEVATORS

*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*



**ANA-003-145**  
Fig. 1L  
Seldin



**ANA-003-146**  
Fig. 1R  
Seldin



**ANA-003-147**  
Fig. 4L  
Seldin



**ANA-003-148**  
Fig. 4R  
Seldin



**ANA-003-149**  
Fig. 51  
Curtis



**ANA-003-150**  
Fig. 52  
Curtis



**ANA-003-151**  
Fig. 11  
Stout



**ANA-003-152**  
Fig. 11A  
Stout



**ANA-003-153**  
Fig. 12

# ROOT ELEVATORS



*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*



**ANA-003-154**  
Fig. A  
Cogswell



**ANA-003-155**  
Fig. 301  
Apexo



**ANA-003-156**  
Fig. 302  
Apexo



**ANA-003-157**  
Fig. 303  
Apexo



**ANA-003-158**  
Fig. 301W  
Seldin



**ANA-003-159**  
Fig. 304W  
Seldin



**ANA-003-160**  
Fig. 1H  
Heidbrink



**ANA-003-161**  
Fig. 2H  
Heidbrink



**ANA-003-162**  
Fig. 40H  
Heidbrink



**ANA-003-163**  
Fig. 41H  
Heidbrink

ROOT ELEVATORS

*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*

ROOT ELEVATORS



**ANA-003-164**  
Fig. 31  
Cryer



**ANA-003-165**  
Fig. 32  
Cryer



**ANA-003-166**  
Fig. 25  
Cryer



**ANA-003-167**  
Fig. 26  
Cryer



**ANA-003-168**  
Fig. 44  
Cryer



**ANA-003-169**  
Fig. 45  
Cryer



**ANA-003-170**  
Fig. 81A



**ANA-003-171**  
Fig. A1C  
Coupland



**ANA-003-172**  
Fig. A2C  
Coupland



**ANA-003-173**  
Fig. A3C  
Coupland

# ROOT ELEVATORS



*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*



**ANA-003-174**  
Fig. 190

**ANA-003-175**  
Fig. 191



**ANA-003-176**  
Fig. 15  
Woodward



**ANA-003-177**  
Fig. 16  
Woodward



**ANA-003-178**  
Fig. 1C  
Coupland



**ANA-003-179**  
Fig. 2C  
Coupland



**ANA-003-180**  
Fig. 3C  
Coupland



**ANA-003-181**  
Fig. 11M  
Mac Millan



**ANA-003-182**  
Fig. 12M  
Mac Millan



**ANA-003-183**  
Fig. 13M  
Mac Millan



**ANA-003-184**  
Fig. 31F  
Friedman



**ANA-003-185**  
Fig. 32F  
Friedman

ROOT ELEVATORS

*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*



**ANA-003-186**  
Fig. 1

3.2 mm



**ANA-003-187**  
# 2  
Kopp



**ANA-003-188**  
# 3  
Kopp



**ANA-003-189**  
# 4  
Kopp

4 mm



**ANA-003-190**  
Fig. 2



**ANA-003-191**  
Fig. 2  
Flohr



**ANA-003-192**  
Fig. 3  
Flohr

3.5 mm



**ANA-003-193**  
Fig. 7  
Cryer



**ANA-003-194**  
Fig. 8  
Cryer

# ROOT ELEVATORS



*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*



**ANA-003-195**  
Fig. 1  
Bein



**ANA-003-196**  
Fig. 2  
Bein



**ANA-003-197**  
Fig. 3  
Bein



**ANA-003-198**  
Fig. 4  
Bein



**ANA-003-199**  
Fig. 5  
Bein



**ANA-003-200**  
Fig. 6  
Bein



**ANA-003-201**  
Berten



**ANA-003-202**  
Berten

ROOT ELEVATORS

*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*



**ANA-003-203**  
Thomson



**ANA-003-204**  
# 1  
Thomson



**ANA-003-205**  
# 2  
Thomson



**ANA-003-206**  
# 3  
Thomson



**ANA-003-207**  
# 4  
Thomson



**ANA-003-208**  
# 5  
Thomson



**ANA-003-209**  
# 6  
Thomson



**ANA-003-210**  
# 1  
Schlemmer



**ANA-003-211**  
# 2  
Schlemmer



**ANA-003-212**  
# 3  
Schlemmer

# ROOT ELEVATORS



*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*



**ANA-003-213**  
# 1



**ANA-003-214**  
# 2



**ANA-003-215**  
# 3



**ANA-003-216**  
# 4



**ANA-003-217**  
# 1



**ANA-003-218**  
# 2



**ANA-003-219**  
# 4



**ANA-003-220**  
# 5

ROOT ELEVATORS

*Елеватори екстракційні*  
*Wurzelheber*

*Élevateurs à racines*  
*Elevadores para raíces*

Atraumatic Extraction Forceps & Root Tip Picks



**ANA-003-221**  
Warwick Root Tip  
Elevator - Left



**ANA-003-222**  
Warwick Root Tip  
Elevator - Right



**ANA-003-223**  
Heidbrink H2-H3  
Double-Ended

Set of Warwick Root Tip Elevator



**ANA-003-224**  
Heidbrink H1  
Single-Ended



**ANA-003-225**  
Heidbrink H2  
Single-Ended



**ANA-003-226**  
Heidbrink H3  
Single-Ended



# ROOT ELEVATORS



*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raices*

Atraumatic Extraction Salvin Elevators & Socket Scoops  
Socket Scoops For Cleaning Out Extraction Sites



**ANA-003-227**  
Misch Socket Scoop  
4mm Angled Left



**ANA-003-228**  
Misch Socket Scoop  
4mm Angled Right



**ANA-003-229**  
Potts Elevator-  
Left



**ANA-003-230**  
Potts Elevator-  
Right



**ANA-003-231**  
Potts "T-Handle"  
Elevator-Left



**ANA-003-232**  
Potts "T-Handle"  
Elevator-Right

*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raices*

Atraumatic Extraction Salvin Elevators



**ANA-003-233**  
DEL 301



**ANA-003-234**  
DEL 302



**ANA-003-235**  
DEL 303



**ANA-003-236**  
DEL 34



**ANA-003-237**  
DEL 34S



**ANA-003-238**  
DEL 36

# ROOT ELEVATORS



*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*

Atraumatic Extraction Salvin Elevators



**ANA-003-239**  
DEL 46R



**ANA-003-240**  
DEL 44



**ANA-003-241**  
DEL 45



**ANA-003-242**  
Cogswell



**ANA-003-243**  
DEL 77R

*Елеватори екстракційні*  
*Wurzelheber*

*Élévateurs à racines*  
*Elevadores para raices*

Atraumatic Extraction Serrated Tip Elevators



**ANA-003-244**  
Serrated Elevator  
2.2mm



**ANA-003-245**  
Serrated Elevator  
3mm



**ANA-003-246**  
Serrated Elevator  
4mm



**ANA-003-247**  
Serrated Elevator  
4mm Left



**ANA-003-248**  
Serrated Elevator  
4mm Right

# ROOT ELEVATORS



*Елеватори екстракційні  
Wurzelheber*

*Elévateurs à racines  
Elevadores para raices*



**ANA-003-249**  
# 1

**ANA-003-250**  
# 2



**ANA-003-251**  
# 3



**ANA-003-252**  
# 1

**ANA-003-253**  
# 2



**ANA-003-254**  
# 1



**ANA-003-255**  
# 3



**ANA-003-256**  
# 4



**ANA-003-257**  
# 2

ROOT ELEVATORS

*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*



Walter F. Barry



**ANA-003-258**  
Fig. 320  
Walter F. Barry



**ANA-003-259**  
Fig. 321  
Walter F. Barry



**ANA-003-260**  
Fig. 322  
Walter F. Barry



**ANA-003-261**  
Fig. 323  
Walter F. Barry



**ANA-003-262**  
Fig. 350  
Walter F. Barry



**ANA-003-263**  
Fig. 351  
Walter F. Barry



**ANA-003-264**  
Fig. 370  
Walter F. Barry



**ANA-003-265**  
Fig. 371  
Walter F. Barry



Lecluse



**ANA-003-266**  
Fig. 1LC  
Lecluse

# ROOT ELEVATORS



*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*



Winter



**ANA-003-267**  
Fig. 11LW  
Winter



**ANA-003-268**  
Fig. 1RW  
Winter



**ANA-003-269**  
Fig. 11LW  
Winter



**ANA-003-270**  
Fig. 11RW  
Winter



**ANA-003-271**  
Fig. 12LW  
Winter



**ANA-003-272**  
Fig. 12RW  
Winter



**ANA-003-273**  
Fig. 14LW  
Winter



**ANA-003-274**  
Fig. 14RW  
Winter

ROOT ELEVATORS

Елеватори екстракційні  
Wurzelheber

Élévateurs à racines  
Elevadores para raíces

Root-tip Picks  
Root-Рада вибірка



**ANA-003-275**  
Fig. 1  
Heidbrink



**ANA-003-276**  
Fig. 2  
Heidbrink



**ANA-003-277**  
Fig. 3  
Heidbrink



**ANA-003-278**  
Fig. 1  
Warwick-James



**ANA-003-279**  
Fig. 2  
Warwick-James



**ANA-003-280**  
Fig. 3  
Warwick-James



**ANA-003-281**  
Fig. 1  
Heidbrink



**ANA-003-282**  
Fig. 2  
Heidbrink



**ANA-003-283**  
Fig. 3  
Heidbrink

## SYNDESMOTOMES



**ANA-003-284**  
Fig. 1  
Chompret



**ANA-003-285**  
Fig. 2  
Chompret



**ANA-003-286**  
Fig. 3  
Chompret



**ANA-003-287**  
Fig. 4  
Chompret

# ROOT ELEVATORS



*Елеватори екстракційні*  
*Wurzelheber*

*Elévateurs à racines*  
*Elevadores para raíces*

## SYNDESMOTOMES



**ANA-003-288**  
11 cm



**ANA-003-289**  
Fig. 0



**ANA-003-290**  
Fig. 1



**ANA-003-291**  
Fig. 2



**ANA-003-292**  
Fig. 3



**ANA-003-293**  
Fig. 4



**ANA-003-294**  
Fig. 5



**ANA-003-295**  
Fig. 6



**ANA-003-296**  
Fig. 7



**ANA-003-297**  
Fig. 8



**ANA-003-298**  
Fig. 9



**ANA-003-299**  
Fig. 11

ROOT ELEVATORS

*Елеватори екстракційні*  
*Wurzelheber*

*Élévateurs à racines*  
*Elevadores para raíces*

## SYNDESMOTOMES



**ANA-003-300**  
Fig. 21



**ANA-003-301**  
Fig. 22



**ANA-003-302**  
Fig. 23



**ANA-003-303**  
Fig. 24



**ANA-003-304**  
Fig. 25



**ANA-003-305**  
Fig. 26



**ANA-003-306**  
Fig. 27



**ANA-003-307**  
Fig. 28



**ANA-003-308**  
Fig. 29

## PERIOTOME HANDLE & BLADES



**ANA-003-309**  
Single end  
6", 15cm



**ANA-003-310**  
P1, 2mm, Straight  
1-1/2", 4cm,



**ANA-003-311**  
P2, 1.5mm, Straight  
1-1/2", 4cm



**ANA-003-312**  
P3, 1.5mm, Angled  
1-3/4", 4.5cm

### PERIOTOME HANDLE & BLADES

- Single end handle can be tapped with a mallet.
- Blades Sold Separately